


Basic Examples

- ▶ cube.txt simple Cube
- ▶ PlotViz3_xml_simple_input.pviz or PlotViz3_simple_input.txt
- ▶ 10Kx10Ksampled.nw.4D.sqrt.sammon.allclusters.pviz 10K Proteomics COG Data

Synchronized Camera View

Synchronized
view



- ▶ Side-by-side analysis for multiple plots
 - Cameras in multiple plots (more than 2) can move in a synchronized way
 - Maintain same rotation angle, camera distance, ...
 - Turn on Features > Camera > Use SyncCamera